







Drill System	Fix 8™ Drill Bit	Ø1.5 Screw (1.2mm drill)	
<u>Standard low speed drills:</u> Aesculap® - Stryker® Anspach® - 3M® Small MicroAire® - Zimmer® (Hall®)	 Standard J-hook Drill Bit	F8-DBS2-04 F8-DBS2-06 F8-DBS2-10	4mm long 6mm long 10mm long
<u>Long handpieces:</u> Aesculap® Hall UltraPower® (Angled Handpiece)	 Long Drill Bit	* DBU2-04 DBU2-06 DBU2-10	4.5mm long 6mm long 10mm long
<u>Hand Drill:</u> Use with small Fix 8™ screwdriver handle	 Hand Driver Drill Bit	F8-DBZ2-04 F8-DBZ2-06	4.5mm long 6 mm long
<u>Hand Drill:</u> Use with large Bioplate® screwdriver handle	 Hand Driver Drill Bit	* DBM2-04 DBM2-06	4.5mm long 6 mm long
Midas Rex® Classic "B" Attachment or "DG2" Attachment	 Use with Midas Rex® "B" Attachment or Order drills from Midas Rex® for "DG2"	* DBB2-04 DBB2-06 DBB2-10	4.5mm long 6mm long 10mm long
<u>Midas Rex III®</u> GS and GS20	 Use with Midas Rex® "S4" Tube	* DBR2-04 DBR2-06 DBR2-10	4.5mm long 6mm long 10mm long

Aesculap® is a trademark of Aesculap, Inc. Anspach® is a trademark of The Anspach Companies. A.O. Synthes® is a trademark of A.O. Synthes, Inc. Bioplate® is a trademark of Bioplate, Inc. Codman® is a trademark of Johnson & Johnson Professional, Inc. MicroAire® is a trademark of MicroAire Surgical Instruments, Inc. Midas Rex® is a trademark of Midas Rex Pneumatic Tool, Inc. Stryker® is a trademark of Stryker Corporation. Zimmer®, Hall® and UltraPower® are trademarks of Zimmer, Inc.

*: Not on contract

IOP, Inc.
 3184 Airway Avenue, Suite B, Costa Mesa, CA 92626
 Tel 800.535.3545 / 714.549.1185 Fax 714.549.0557
 www.iopinc.com

FIX 8 is a trademark of IOP, Inc. Manufactured by Bioplate®, Inc. ©2005

Drill System	Fix 8™ Drill Bit	Ø1.5 Screw (1.2mm drill)	
<u>Anspach Black Max®</u>	 Black Max Type	*	DBT2-04 4.5mm long DBT2-06 6mm long DBT2-10 10mm long
<u>Aesculap</u>	 Aesculap High Speed Type	*	DBA2-04 4.5mm long DBA2-06 6mm long DBA2-10 10mm long
<u>Anspach Micro Max®</u>	 Micro Max Type	*	DBF2-04 4.5mm long DBF2-06 6mm long
<u>Osteomed</u>	 Osteomed Type	*	DBG2-04 4.5mm long DBG2-06 6mm long DBG2-10 10mm long
<u>Stryker TPS</u>	 Stryker TPS Type	*	DBQ2-04 4.5mm long DBQ2-06 6mm long
<u>Hall UltraPower®</u> (High Speed Handpiece)	 Use with Hall UltraPower® High Speed Handpiece Available from Bioplate®	*	DBH2-04 4.5mm long DBH2-06 6mm long DBH2-10 10mm long
<u>Quick Coupling:</u> AO Synthes® MicroAire® Contra Angle Stryker® Contra Angle	 Mini Quick Coupling Drill Bit Available from Bioplate®	*	DBD2-04 4mm long DBD2-06 6mm long DBD2-10 10mm long

*: Not on contract