

SutureGroove

Sterile Gold Eyelid Weight Implants **slim™**

New slim-style gold eyelid weights have arrived.

SutureGroove™ now offers a complete selection of eyelid weight designs for the surgical treatment of lagophthalmos.

● Improved Cosmesis & Patient Satisfaction

SutureGroove Slim™ weights are 40% narrower than the traditional design, providing a more-contoured, less conspicuous appearance within the eyelid.

● Innovation & Design



The SutureGroove™ design features recessed suture channels that allows for a more stable and secure fixation of the suture. This channel embeds the suture, achieving a stable, virtually immovable position.

● Sterile & Secure Packaging

Each weight is packaged sterile with easy “Peel and Stick” tracking labels for effortless patient recording.

● Manufactured Precision

All SutureGroove™ weights are forged from pure gold. Each weight is quality-inspected for a smooth surface and rounded, tapered edges. Each implant is weight-verified for accuracy.



SutureGroove Slim Weights

Item #	Weight
GSLM-08	0.8 gram
GSLM-10	1.0 gram
GSLM-12	1.2 gram
GSLM-14	1.4 gram
GSLM-16	1.6 gram



3151 Airway Avenue, Suite 1-1
Costa Mesa, CA 92626
(714) 549-1185
Fax: (714) 549-0557

Customer Service:
(800) 535-3545

www.iopinc.com

© 2002 IOP, Inc.
SutureGroove Slim is a trademark of AMI is used with permission.